

Replacement Dental Floss 2 pcs. Yaweco 2pcs x 25

EAN: 4280001065277 FABRICANTE: YAWECO



Dental floss replacement. Natural silk dental floss.

DESCRIPTION

Replacement of 2 rolls of natural dental floss for the Yaweco dental floss case (reference Y110010).

Dental floss with 100% natural silk, without synthetic ingredients.

Helps to deeply clean the interdental space, removing plaque without damaging the gums.

Biodegradable thread and packaging. 25 metres each refill.

About the composition of the Natural Dental Floss Refill 2 Rolls

- Natural silk: soft and highly effective material thanks to its porosity.
- Coated with beeswax: makes it easier to apply dental floss between the mouthpieces. Excellent interdental cleaning.

How to use Natural Dental Floss Refill 2 Rolls?

How to refill the dispenser:

Once you have finished flossing, open the dispenser by separating the top part from the bottom. The two parts fit together but can be easily separated by applying a little pressure.

Remove a roll from the refill pack.

At the bottom of the dispenser is the spool of thread that has run out. Remove it, the new roll has its own spool.

Insert the new spool and guide the line as indicated on the replacement pack. Pass the end of the line through the notch to the left of the blade and around it. Close the dispenser by pressing the top onto the bottom.

Prices valid except for changes on our website. (IVA/VAT) included. Doc. Version bf-trf-G-20241210
<https://bio-farma.es/buy/f/yaweco/replacement-dental-floss-2-pcs-yaweco-2pcs-x-25>





Replacement Dental Floss 2 pcs. Yaweco 2pcs x 25

EAN: 4280001065277 **FABRICANTE:** YAWECO



Make sure the two flat surfaces on the top and bottom match so you can rotate the dispenser to open or close the thread outlet and blade.

If the two surfaces do not match exactly, reopen the dispenser, readjust it, and close it again.

Prices valid except for changes on our website. (IVA/VAT) included. Doc. Version bf-trf-G-20241210
<https://bio-farma.es/buy/f/yaweco/replacement-dental-floss-2-pcs-yaweco-2pcs-x-25>

<https://bio-farma.es> - info@bio-farma.es

C/ Vermell, 20 3r 3a - 17252
St. Antoni de Calonge (Girona) Spain

 +34 696 79 77 77

